




CPAPMAN.COM
BEST SERVICE SINCE 1997

Click Here To Order Today
Sizing guide

Instructions:


- 1 Print out this Adobe PDF File at 100% on your printer. Use the print accuracy check at right to ensure printing at the proper size.
- 2 Fold the bottom part of this page on the dashed line.
- 3 Cut along the dotted line around the fitting gauge.
- 4 Place the fitting gauge on your nose as shown in the picture below.
- 5 Determine the appropriate size by which area your nose reaches.
(S – Small, M – Medium or L - Large)

Philips Respironics always recommends consulting a physician or registered Sleep Therapist in fitting a mask to ensure a proper sizing and fit.



FOLD along dashed line

Print accuracy check



Fitting gauge

CUT along dotted line

