

Pilot Flex - Run Times

Brand	Machine Model	Pressure (cm H ₂ O)	Humidifier and Heated Tubing	Run Time
ResMed	AirMini	Medium - 10	No	14 hrs 16 min
	AirSense 10		Yes	2 hrs 16 min
			No	11 hrs 36 min
	AirSense 11		Yes	2 hrs 17 min
			No	12 hrs 48 min
Philips	DreamStation 1	Medium - 10	Yes	3 hrs 2 min
			No	18 hrs 47 min
	DreamStation 2		Yes	2 hrs 53 min
			No	19 hrs 2 min
Löwenstein	Löwenstein Prisma SMART	Medium - 10	No	12 hrs 1 min
BMC	Luna G3	Medium - 10	Yes	3 hrs 6 min
			No	14 hrs 47 min